PROVI MASTER F

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PRODUCT DESCRIPTION:

PROVIMASTER F is a radiopaque light-cured resin for temporary restorations of dental preparations, inlays, onlays, crowns, endodontic treatments, and sealing of access holes in implant-supported prostheses. It is intended for procedures that require tooth protection before the final restoration, keeping it properly prepared. The components of PROVIMASTER F provide adequate strength to temporarily withstand masticatory wear. Supplied in a single syringe, it is easy to handle and apply. Once polymerized, the resin allows for easy removal due to its elasticity.

INDICATIONS: PROVIMASTER F is indicated for temporary filling techniques, sealing of endodontic access, and use on implant components, preserving the tooth's shape and structure.

PRESENTATION: Available in individual packaging: 1 syringe of 3g in the following shades – Universal, A1, and A3 – as indicated on the box along with the instructions for use.

Also available in a kit containing 1 syringe of each shade. (Wilcos follows the VITA Classical Shade Guide standard.)

COMPOSITION:

UDMA, TEGDMA, EDAB, BHT, sodium fluoride, photo initiator, and fillers.

APPLICATION TECHNIQUE: Begin with proper cleaning (asepsis) and drying of the cavity. Apply a single layer of PROVIMASTER F into the prepared cavity. Shape the resin using a suitable dental

instrument. To preserve tooth mobility, keep one contact point free of material. Light-cure for 20–40 seconds using a halogen light at approximately 1200 mW/cm². If the curing tip cannot be placed close to the resin, increase the exposure time by another 20–40 seconds. If needed, use fine silicone tips for polishing.

WARNING, PRECAUTION, RESTRICTION:

WARNING: Use a curing light with at least 1200 mW/cm². Lower-power units may require longer curing time. This product is intended for professional use only in dental practice, and must follow the manufacturer's instructions. The resin may adhere to other dental materials or cavity liners; isolation with glycerin or petroleum jelly is recommended. Avoid materials containing eugenol or thymol, as they inhibit curing. No significant color change is expected during use, but regular dental monitoring is recommended. All components are biocompatible with oral tissues and fluids. Do not use expired product or in patients with known hypersensitivity to any of its components.

PRECAUTION: Ensure proper asepsis and curing time with the appropriate light intensity. Monitor for any patient discomfort due to heat during curing. Make sure no excess material remains in the pulp cavity. For temporary restorations lasting over 2 weeks or placed on occlusal/marginal surfaces, reinforce with composite resin. Apply a small amount to the desired area, shape, and cure with the resin. The material slightly expands after curing, promoting excellent marginal adaptation. In cases of multiple adjacent restorations, remove excess material before light-curing.

RESTRICTION: Ingestion may cause gastric irritation. Single-use product – once cured, it cannot be reused or reverted to its original state. Do not reprocess. Avoid exposure to temperatures above 25°C or direct sunlight.

CONTRAINDICATIONS: Ingestion may cause gastric irritation. Do not reuse the product once polymerized. Do not reprocess. Avoid exposure to temperatures above 25°C or direct sunlight.

ADVERSE EFFECTS: No adverse effects related to the product or its components have been reported to date.

EXPIRATION DATE: 36 months, provided it remains in original, sealed packaging and stored according to instructions. Do not use after the expiration date.

STORAGE: Store at room temperature (15°C to 30°C) in a dry place, protected from light and heat, and in its original packaging. Light exposure may alter the formulation and affect product properties.

DISPOSAL: PROVIMASTER F is for single use only and must be discarded according to local regulations.

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